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Revision 04

**USER ALERT**

- Before using the product, verify that the version of this digital Instructions for Use matches the version indicated on the labeling of the acquired product.

- A printed version of this Instructions for Use is available free of charge, including shipping. To request a printed copy, please contact our Customer Service at [sac@medicone.com.br](mailto:sac@medicone.com.br) or call 0800 722 2728.

**INTENDED USE**

The Medicone Testicular Implant was developed as a cosmetic solution or corrective device for patients who have experienced loss, malformation, absence, or disease affecting the natural testicle.

**PRODUCT DESCRIPTION AND COMPOSITION**

The Medicone Testicular Implant consists of a medical-grade silicone shell that is chemically and mechanically resistant. Volumetric comfort is ensured by a silicone filling similar to human tissue, shaped in a mold to match the form, density, and overall consistency of natural tissue.

**COMMERCIAL PRESENTATION**

The Medicone Testicular Implant is supplied individually and packaged in a dual triple-barrier sterile packaging system.

**PRODUCT INDICATIONS**

- General indications for the use of a testicular implant;
- Patients who had a natural testicle removed due to disease, trauma, or accident;
- Children with congenital malformation or absence of one or both natural testicles;
- Patients with anorchia, cryptorchidism, testicular atrophy, or other genetic or infectious causes.

**INSTRUCTIONS FOR USE**

- Check if the implant shows any signs of leakage.
- If any damage is detected, the product must not be used.
- Prepare the patient and surgical environment according to routine guidelines. Proceed with anesthesia, vital sign monitoring, and proper patient positioning.
- The surgeon should select the surgical approach deemed most appropriate for positioning the implant.
- To secure the implant, the surgeon must use the fixation site located at the top of the implant.

**WARNING:** Do not puncture or suture the implant.

**STORAGE CONDITIONS**

- Keep the packaging sealed and intact until use.
- Store at room temperature, protected from heat and humidity.

**STERILIZATION**

The Medicone Testicular Implant is supplied sterile by ethylene oxide (ETO). The sterilization validity period is three to five years from the sterilization date, which is clearly indicated on the product packaging.

**SINGLE-USE PRODUCT – REUSE PROHIBITED**

The Medicone Testicular Implant is intended for single use only. Reuse or re-sterilization is not recommended, as it may compromise product performance and safety, and pose a risk of infection.

## IMPLANT SHELF LIFE IN THE BODY

Medicone defines an expected implantable life of ten years. However, this period may be adjusted based on clinical justifications. Annual medical follow-up is recommended for patients using the Medicone Testicular Implant.

## PRECAUTIONS

- The implant surface must not come into contact with talcum powder, ointments, or oils. Sterile surgical gloves should be rinsed with sterile saline solution before handling the implant.
- The implant is supplied sterile. To maintain its sterility, handling must be minimized. Use of antiseptics, disinfectants, or cleaning agents on the implant is strictly prohibited.
- The selection and use of this product must be performed by trained medical professionals who can determine the most appropriate size and type of implant for the patient.
- The physician must assess the patient's health status and evaluate the feasibility of the procedure.
- Hospital medical device. The manufacturer recommends single use only. Reuse may lead to patient infection.
- The package must be intact, and the product within its expiration date, before use.

## POSSIBLE COMPLICATIONS

- Possible complications associated with implant use should be discussed with the patient prior to surgery.
- Complications may include risks related to medications, surgical techniques, intolerance or allergic reactions.
- Signs of infection may include pain, heat, and swelling at the implant site. In such cases, patients may not respond to treatment and implant removal may be required.
- Fluid accumulation at the implant site is common postoperatively.
- In rare cases, implant migration or poor positioning may be linked to the surgical technique used.
- Formation of a fibrous capsule around the implant is a normal physiological response after surgery.

This product must be indicated and handled only by trained and qualified professionals. Any questions or adverse events should be reported to Medicone. Medicone is not liable for accidents or consequences resulting from improper use of the product.

Although the implant consists of an inert shell, it has high elasticity and may attract small external particles such as metal shavings. Handling under contamination risk conditions should be avoided.

## ADVERSE EFFECTS

It is the physician's responsibility to inform the patient that a future surgical procedure may be required to remove or replace the implant. The following adverse events have been reported in the medical literature:

- **Calcification** – May cause discomfort; observed in areas adjacent to the implant.
- **Capsular Contracture** – Of unknown cause; believed to be reduced with proper surgical technique and careful handling. May occur when the implant is placed subcutaneously, exerting pressure on the surface.
- **Implant Displacement** – May result from patient movement or desired positioning; more common in posterior fixation techniques.
- **Pain** – Some patients may experience pain or hypersensitivity, which should always be investigated. Inadequate surgical techniques may intensify pain.
- **Extrusion** – May result from excessive pressure or trauma at the implant site and/or interrupted healing.
- **Hematoma** – May lead to infection or seroma. Should be resolved or removed. Meticulous hemostasis is the key preventive measure.
- **Infection** – Often due to poor aseptic technique. Typically confirmed by local culture and treated with antibiotics, drainage, and possibly implant removal.
- **Skin Necrosis** – May occur when implants exceed the anatomical space, leading to skin tension and necrosis.
- **Loss of Sensation** – Usually temporary and reversible over time but may persist for 2 to 3 months.
- **Immunological Response** – Some cases reported fibrous or hardened tissue due to implant presence, but no evidence of immunological rejection has been observed.

Proper indication and use of the implant must be performed by specialized medical teams. If required, the implant should be replaced with one of appropriate size and characteristics, as defined by technical specifications.

**IMPLANT RUPTURE OR DAMAGE**

During surgery, the implant may be damaged by sharp surgical instruments. If rupture or perforation is suspected, the implant should be carefully evaluated. A clinical decision must then be made regarding corrective actions, including immediate replacement, if necessary.

**DISPOSAL**

After use, all components must be discarded in accordance with institutional protocols and applicable biosafety regulations for potentially contaminated materials.

**STERILE PRODUCT – Sterilized by Ethylene Oxide (ETO)**

**SINGLE-USE PRODUCT – DO NOT REPROCESS**

**Manufacture Date, Expiration Date, and Lot Number:** See product packaging.

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